

# **PennTech Newsletter**

NL# 716

Date: November 2, 2007

## Quote of the Week:

*“Begin somewhere; you cannot build a reputation on what you intend to do”*

...Liz Smith

## Next Week's Topic:

- Proper Vial Cooling After Sterilization -

PennTech Machinery Corp.  
103 Steamwhistle Drive  
Ivyland, PA 18974  
USA  
Phone: 215.396.2200  
Fax: 215.396.6774  
info@penntech-corp.com




Author: Ger Smit  
gsmit@penntech-corp.com

## **PennTech Complete Aseptic Production Lines**

During the next Interphex Show (March 26-28, 2008 in Philadelphia) we are planning to introduce a complete PennTech aseptic production line. Vial washing - sterilizing/ depyrogenation - filling, stoppering, capping – trayloading – completely integrated.

We are now able to offer complete, compact, high quality, aseptic production lines.....and the US-dollar is helping. Our prices are 30-50% lower than our European colleagues. Look at our documentation: Complete VMP (Validation Master Program) to help the customer in expediting the validation process, 3-D exploded assembly drawings of the equipment to simplify maintenance and parts ordering, material certifications, welding certifications. PennTech has received rave reviews from its customers concerning the quality and comprehensiveness of its documentation.



### **Specifications:**

Vial Range 1-100ml

Throughput: up to 6000vials/hr (2ml vials)

Vial Washing: model RW-500

Sterilizing/Depyrogenation: model PST-450/2340

Filling/Stoppering/Capping: model FSC-222

Trayloading: model TL-100

***PennTech...uniquely designed to outperform***